

Aqualog Version 3.6.10.1: Release Notes

Product Release Information

Product	Aqualog Software
Release Number	3.6.10.1
Release Date	12/22/2014
Customer Support	For more information or support, please visit our website

Introduction

This is a release to add features to version 3.6.8. Your feedback is important to us. Please help us make this the best product possible. Keep in mind that we are continuing to work on Aqualog and things may change in the future.

What's New

1. Continuous Sample Queue blank file subtract. This feature allows the user to select a previously collected (known good) blank file.
2. Support for 2D graph overlay tool.
3. Fixes for memory leak issues.
4. Fix for QSU Normalization tool.

Installation and Upgrade Notes

Patch can be applied to Aqualog 3.6.3 or newer versions.

Recent Changes

Type	Tracking #	Description
feature		Continuous Sample Queue blank file subtract. This feature allows the user to select a previously collected (known good) blank file.
feature		Support for 2D graph overlay tool.
Fix		Fixes for some memory leak issues.
Fix		Fix for QSU Normalization tool.

Known Problems and Workarounds

Type	Tracking #	Description	Workaround
issue		<i>Sample Queue</i> may require user to Create new default experiment.	User must create new default experiment.

